


EXPERIENCE VISEON



MaxView™ System – Posterior Solution

Maximizing Minimally Invasive Spine Surgery

Viseon MaxView™, the only high-definition, single-use surgical camera for minimally invasive spine surgery, enables intraoperative visualization of the surgical field on a head's up display, in one easy-to-use system.

The MaxView™ Posterior Solution is designed for use in posterior approaches to the spine and provides:

- An unobstructed view down retractor tube for better visualization of anatomy
- Head's up display on operating room monitors for optimized ergonomics and OR workflow
- Enhanced teaching experience to improve learning for residents and fellows

 **VISEON**
INC.

MaxView™ System – Posterior Solution

Maximizing Minimally Invasive Spine Surgery

The **MaxView™ Posterior Solution** minimizes disruption and maximizes efficiencies in minimally invasive spine surgery while providing exceptional visualization. Features include:

- Excellent image quality and highly magnified views of the surgical field
- Retractor tube with integrated LED lighting to eliminate the need for a separate light source
- Slim profile design that is easy to set up and has minimal operating room footprint
- Imager and tube that are sterile, single use to ensure a sharp image and sterility for every case
- Universal connection to operating room monitors and compatible with most table arms and dilators

MaxView™ System – Posterior Solution includes:

MODEL NUMBER	DESCRIPTION
Refer to Tube size chart below	Posterior Set , includes sterile, single use Viseon Tube and Viseon Imager
ICB-02	Accessories , non-sterile, reusable Image Control Box (ICB) and connecting cables

Posterior Set, Tube size chart*:

TUBE SIZE	OUTER DIAMETER	
	18mm	22mm
LENGTH		
40mm	VSY-1840	VSY-2240
50mm	VSY-1850	VSY-2250
60mm	VSY-1860	VSY-2260
70mm	VSY-1870	VSY-2270
80mm	VSY-1880	VSY-2280
Dilator Compatibility	16.8mm	20.8mm

*The Viseon Tubes are available in multiple diameters and lengths to accommodate a variety of patient anatomies and spinal surgical procedures.



Posterior Set

Image Control Box

For additional product information or to request a demo, please contact Customer Service.

Please refer to the instructions for use for a complete list of contraindications, warnings, and precautions.

Caution: Federal law restricts this device to sale by or on the order of a physician.

